

KF Application Note No. K- 3

Title:	Water in ammonium and potassium peroxodisulphate (persulphates)
---------------	--

Summary:	The water content of ammonium and potassium peroxodisulphate is determined according to Karl Fischer using two-component reagents. To prevent unwanted side reactions the determinations are carried out at -20 °C. Because the potassium salt is insoluble in the solvent, a high-frequency mixer is used to disintegrate the salt particles.
-----------------	--

Sample:	Ammonium and potassium peroxodisulphate
----------------	---

Sample Preparation:	none
----------------------------	------

Instruments and Accessories:	701 KF Titrino or 720 KFS Titrino, 703 Titration Stand, printer, Polytron PT 1200 Disintegrator, low temperature circulation system
-------------------------------------	---

Analysis:	Wait for a steady drift below 6 uL/min, then add ca. 1 ... 1.5 g sample using a glass weighing spoon. An extraction time of 2 min with intensive stirring (Polytron PT 1200, speed «3») has been used for the automatic determination.
------------------	---

Reagents:

Solvent: 50% formamide + 50% Hydranal Solvent (Riedel-de Haën)

Titant: Hydranal Titant 5 (Riedel-de Haën)

Results:	(NH ₄) ₂ S ₂ O ₈ : AVG(5) = 377 +/- 41 ppm water K ₂ S ₂ O ₈ : AVG(5) = 419 +/- 28 ppm water
-----------------	---

Settings:	701 KF Titrino
	>titration parameters
	extr.time 120 s
	stop crit.: drift
	stop drift 20 uL/min
	>preselections
	conditioning: on
	req.smpl size: on
	report: full